

FASTER

FASTER - TIMI 24

Protocol Design



ST Elevation MI < 6 h

ASA

TNK-tPA 0.53 mg/kg

Tirofiban bolus + infusion
(10 to 15 µg/kg) 0.15 µg/kg/min

STD Heparin

↓ TNK-tPA dose

Low Dose Heparin

Primary Endpoint Part A (Dose Finding): **TIMI 3 Flow at 60 min**

Primary Endpoint Part B (Dose Confirmation): **ST segment resolution at 60 min**


FASTER

FASTER - TIMI 24


Results



TIMI Flow at 60 min



30 d Death, MI, Urgent Revasc



*Corrected for multiple events